

11.12.19

Objectives

- Seed 1.75ul plates using OP50-GFP
- Seed more 6cm plates
- Pick AML-283/284 and N2
- Check plate stocks
- Pour more plates

Observations

- At about 41hrs plates AML-283/284 and N2 all seem to be at the same stage - indicates that divergence begins with vulval development but will follow up on it tomorrow
- OP50-GFP plate did not form many single colonies because plate only dried for 10 minutes before streaking - allow to dry for 1hr this time before re-streaking

From:

<http://wiki.scholz-lab.com/> - Scholz Lab Wiki

Permanent link:

<http://wiki.scholz-lab.com/doku.php?id=wiki:documentation:wiki:documentation:nicolina:11.12.19>

Last update: **2019/12/12 10:10**

